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OXYGEN

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SAFE USE IN MEDICAL FACILITIES

Kim Trinh Nguyen – Secretary General

Asia Industrial Gases Association in Vietnam

17th Floor, Prime Centre building, 53 Quang Trung str., Hai Ba Trung Dist., Ha Noi, Vietnam

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DISCLAIMER

All publications of AIGA VN or bearing AIGA VN's name contain information, including work instruction, safety procedures applying, technical information, legal requirement, events that were obtained from sources believed by AIGA VN to be reliable and/ or based on technical, safety information and experience currently available from members of AIGA VN at the date of issuing.

AIGA VN do not ensure nor accept any liability as to the accuracy, completeness or correctness of the information contained in these publications. While AIGA VN recommends that the user refer to or apply it into their working, such reference to or apply thereof by the user is purely voluntary and not binding.

AIGA VN has no control whatsoever as regards, performance or nonperformance, misinterpretation, proper or improper use of any information or suggestions contained in AIGA VN's publications by any person or entity (including AIGA VN members). However, AIGA VN recommends that users do not modify, change the content and use publications of AIGA VN for improper purposes, and AIGA VN expressly disclaims any related liability.

AIGA VN's publications are subject to periodic review and the users are cautioned to obtain the latest edition.

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

Áp suất 1atm

Pressure

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
PERCENTAGE OF OXYGEN IN THE ATMOSPHERE

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
Never use flammable and Electric discharge materials in Oxygen gas equipment.

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Avoid heat, flames, electricity and other sources of ignition

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Such as grease, oil, asphalt, or Gaskets made of flammable materials

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Oxygen is NOT a flammable gas but is a strong oxidizing agent, which is an element that supports on fire and makes a fire more intense

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ASU (Air Separation Unit) produces Nitrogen, Oxygen and Argon with high capacity at factories

Oxygen from ASU is produced by compressing air and gradually cooling it down to the liquefaction temperature of each substances presenting in the air: oxygen (-183°C) Argon (-186°C), and nitrogen (-196°C); and therefore, each substance can be separated


Oxygen product obtained from ASU is cryogenic liquid and has a purity of 99.6%, it then will be filled into the storage tank

Supply to medical facilities using high volumes of oxygen by tanks

Supply to medical facilities using a small volume of oxygen by high pressure cylinders (liquefied oxygen was vaporized into gaseous phase before filling into the cylinders)

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PSA (Pressure Swing Adsorption) Oxygen 93% can be manufactured on-site within medical facilities, under the responsibility of medical facilities qualified personnel

Suitable for production and supply of oxygen at sites where access for cylinders and liquid oxygen supplies are difficult

93% O₂ manufactured via PSA process are allowed for medical use by most country Authorities

Oxygen 93% are produced using pressure swing adsorption process (PSA)

Oxygen product 93% is gaseous and has a purity of 93%, which is filled directly into the cylinder or supplied direct to hospital bed

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O₂ | PSA Oxygen Generator
High Purity Oxygen Production

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THE DIFFERENCE BETWEEN THE 2 PRODUCTION METHODS

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Oxygen produced by ASU – Air Separation Unit is used

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In the industry called industrial oxygen, and used such as welding, metallurgy, rocket launch, sewage treatment, etc.,

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Oxygen produced from PSA with a purity of 93% may only be suitable for patients with mild COVID.

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<https://duocdienvietnam.com/oxygen/>

According to Pharmacopoeia V, the minimum O₂ purity should be 99.5%. On the other hand, many medical facility processes and equipment do not use 93% O₂ but are required to be O₂ 99.6%

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CONVERSION CYLINDERS OF OTHER GASES TO MEDICAL GAS CYLINDER

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When oxygen requirement increases in medical facilities due to the Covid-19 epidemic, the number of medical oxygen cylinders is not enough, it is possible to convert other gas cylinders to medical oxygen cylinders.

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
Liquefaction of oxygen from liquid to gas → filling the cylinder

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HAPPENED INCIDENT AND LESSONS LEARNED

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On 29 Sep 2008, the fire occurred at 08:30 hrs in a four beds intensive care ward in the Cardiac Wing

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Paediatric patients were being treated with Berlin Heart machines (where a lot of oxygen is used)

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The fire was all over in 10 minutes and fortunately no patient was hurt

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
According to AFP news agency on April 21, at a Covid-19 hospital, 22 patients died when the oxygen supply was suddenly interrupted.

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The liquid oxygen tank has a problem, causing a supplying oxygen interrupted about 30 minutes

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At the site, a tank lorry is charging to the storage tank, so the causes may include:

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3

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PIPELINE OXYGEN CONCENTRATION FELL TO 40%

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On 28 Oct 2005, an anaesthetist noticed that the percentage of oxygen in the hospital supply system had dropped from the normal 100% to around 40%, raising independent alarms throughout the hospital.

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Anaesthesiologist mandated no new procedure to commence, and all operating theatres were suspended

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HAPPENED INCIDENT AND LESSONS LEARNED

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Med O2 line was contaminated with Med Air by contractor in the absence of installation drawing

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*Existing lines were **wrongly labelled***

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Contractor didn't check the line content before modifying pipework

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*Hospital Authority was **not notified about commissioning** of new lines*

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A hospital wanted to modify its existing gas supply for an expansion project. The modification was to move existing piping to facilitate civil works (by another contractor) and to add connections to supply gases to 8 delivery rooms and an operating theatre

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When performing pre-startup checks in the operating theatres on the next working day, the anesthetist's stations started to alarm indicating low oxygen level in supplied gas. Fortunately, No patients were involved!

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Mandatory supervision and oversight must be provided for such as work on hospital installations.

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Contractors should be trained and qualified as employees to be able to work unsupervised

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ACCIDENTS HAPPENED IN COVID-19 TREATMENT HOSPITALS IN IRAQ

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According to CNN dated Apr 26, 2021, fire occurred at a hospital in Baghdad

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According to Reuters dated Jul 13, 2021, fire occurred at a hospital Hussain - Iraq

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96% of accidents happen due to human factors

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4% is due to unsafe conditions (nature)

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PERSON WHO WORK WITH OXYGEN SHOULD HAVE

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